

ORP Buffer Solution HI7021L



????????????? ?????????????? ?????????????????? > ORP Buffer Solution HI7021L HANNA ORP Test Solution at 240 mV (@20°C)

?????????: ??????????????????????????
???? :

[??](#)

?????????????: [HANNA](#)

???????????

ORP Buffer Solution HI7021L

Buffer Solution , ORP buffer solution

ORP Solution for Platinum and Gold Electrodes. Reads 240 mV @ 25°C/77°F ORP standard solutions allow you to test the precision of ORP electrodes.

If the reading is outside the indicated interval, clean and condition your ORP electrode in HANNA® pretreatment solution. Use HI 7092 for oxidizing or HI 7091 for reducing pretreatment.

Description	ORP Test Solution at 240 mV (@20°C)
Bottle	500 mL

{tortags,1031,3}

?????????
??.

