#### Accessories

Below is a small selection of the accessories available for the Gas-Pro.



#### 1m reactive gas tubing AC0301

1m reactive gas tubing



#### Gas-Pro CSE Kit Case AC0610

Hard shell case including foam insert with space for Q-Test, 34 or 58ltr cylinder, power and comms cable, drop line, Gas-Pro, flow plate, filters, chest harness and straps and instructions - contents not included.



#### Gas-Pro Hard shell case AC0600

Includes custom foam insert with space for standard items - contents not included.



### Sampling tube



Compatible sampling tube for Gas-Pro. Select length from list and whether standard or low adsorption material for reactive gases. All include tube insert to enable attachment to flow plate.



### Sampling Tube Ball Float Probe AC0512

Connects to sampling tube to protect from water being drawn up when sampling in potentially problematic locations.



#### Tube insert pack AC0500

Replacement inserts to allow connection of sampling tube to Gas-Pro flow plate. Box of 10.



### Elbow tube insert pack AC0511

90 degree inserts to reduce protrusion when connecting sampling tube to Gas-Pro flow plate. Box of 10.



## Hand aspirator bulb AC0504

For manual gas sampling. Use with appropriate flow plate and sampling tube.



### Water trap and spare filters AC0502

Protects detector sensors and pump from particles and moisture. Includes place of 10 replacement filter elements.



## Chest harness plate AC0506

Enables detector to be worn on the chest or over the shoulder. Use with straps AC0507 or AC0508.



### Chest harness straps AC0507

Set of 2 straps which in combination with the chest harness plate enable the detector to be worn on the chest.



### Shoulder Single strap AC0508

Single straps which in combination with the chest harness plate enable the detector to be worn over the shoulder.



### 6m Drop line AC0509

6m Drop line to enable pre-entry checks. Includes clip to attache to D-loop on detector.



# Gas-Pro un-pumped flow plate AC0101

Replacement flow plate for un-pumped Gas-Pro units.



## Gas-Pro Pumped flow plate AC0100

Replacement flow plate for pumped Gas-Pro units.



# Multiregion power lead and supply CH0100



# Gas-Pro charger cradle CH0105

Desk or wall mounted cradle for charging Gas-Pro.



## USB communications lead CH0103

Enables communication with PC (non charging).



# Vehicle charging adapter CH0106

Provides power from 12/24V vehicle outlet sockets. Use with CH0106.